

SZ680 Trinocular Zoom Stereo Microscope with Digital Camera

Product Description

Stereo microscope is also named dissecting microscope. Stereoscopic microscope can provide stereo character with strong visual effect from different angles in one time. Stereo microscope has unparalleled image quality and super wide visual effect.

It can be widely used in the training and education of natural science in universities and research institutes, Daily inspection of medical institutions, bioengineering and scientific research, industrial assembly, testing, and quality control, especially in the IT industry.

The SZ680 trinocular zoom stereo microscope with digital camera can connect a computer, and the user can deal with the images thru software.

Product Feature

- High performance, high quality, and high cost performance
- Adopts dustproof, oil-proof, waterproof, anti-mildew
- Daily inspection of medical institutions, bioengineering and scientific research, industrial assembly, testing, and quality control, especially in the IT industry
- It magnifies the image of object and displays it on the screen of the computer. It can save, enlarge and print the image.
- With measurement software, you can measure various data.



SZ680TP



SZ680T2L

Product Parameter

Product Name	Trinocular Zoom Stereo Microscope	
Model	SZ680TP	SZ680T2L
Objectives Zoom Range	0.68X-4.7X	
Zoom Ratio	1:6.8	
Working Distance	110mm	
Stereo Angle	12°	
Viewing Angle	35°	
Interpupillary Distance	50-75mm (adjustable)	
Video Output Distribution	Standard 0:100, optional 50:50	
Eyepiece	WF10X/23mm, high eyepoint	
Objective	1X	
Focusing Unit	Yes	Yes
Illumination	No	Reflected & transmitting 3W LED
Plate	White & Black plate	White & Black plate Glass plate
Base Stage	330*300mm	330*300mm
C-mount	0.5X C-mount	0.5X C-mount
Digital Camera	USB2.0 3mp CMOS camera	
Optional Accessories	WF15X/16mm Eyepiece	
	WF20X/12mm Eyepiece	
	0.5X additional objective	
	2.0X additional objective	
	1X C-mount	
	0.7X C-mount	
	10X Eyepiece Reticule	
	USB2.0 5MP CMOS camera	
	USB3.0 12MP CMOS camera	



3mp camera



5mp camera



12mp camera

Product Application



About Us



Chongqing Scope Instrument Co., Ltd. is a professional supplier who focus on laboratory instrument, and our team has more than 10 years rich experience in the field of laboratory instrument. We have perfect relationship with most of leading manufacturers in the field of laboratory instrument in China, and supply perfect service and different kinds of quality laboratory instrument from different leading manufacturers to our customers. We now serving for more than 80 customers from Europe, Asia, Africa, South America and Australia. Also, we have established good relationship with more than 200 customers in worldwide market. Our laboratory instrument are widely used in the field of labs, education, medicine, petroleum, geological mineral and chemical industry, and electronic industry.

Product Application



Payment

Western Union

T/T



Shipment



FAQ

Are you a trading company??

Yes, we are a trading company now.

What type of payment do you have available??

We accept TT, Western Union.

Do you accept OEM term?

Yes, we accept OEM.

How is the package?

We pack our goods with foam and hard carton.

